

# इंटरनेट

# मानक

## Disclosure to Promote the Right To Information

Whereas the Parliament of India has set out to provide a practical regime of right to information for citizens to secure access to information under the control of public authorities, in order to promote transparency and accountability in the working of every public authority, and whereas the attached publication of the Bureau of Indian Standards is of particular interest to the public, particularly disadvantaged communities and those engaged in the pursuit of education and knowledge, the attached public safety standard is made available to promote the timely dissemination of this information in an accurate manner to the public.

“जानने का अधिकार, जीने का अधिकार”

Mazdoor Kisan Shakti Sangathan

“The Right to Information, The Right to Live”

“पुराने को छोड़ नये के तरफ”

Jawaharlal Nehru

“Step Out From the Old to the New”

IS 6582 (1971): Bio-assay methods for evaluating acute toxicity of industrial effluents and wastewaters [CHD 32: Environmental Protection and Waste Management]



“ज्ञान से एक नये भारत का निर्माण”

Satyanarayan Gangaram Pitroda

“Invent a New India Using Knowledge”



“ज्ञान एक ऐसा खजाना है जो कभी चुराया नहीं जा सकता है”

Bhartrhari—Nitiśatakam

“Knowledge is such a treasure which cannot be stolen”



BLANK PAGE



IS: 6582 - 1971

( Reaffirmed 1987 )

## *Indian Standard*

# BIO-ASSAY METHODS FOR EVALUATING ACUTE TOXICITY OF INDUSTRIAL EFFLUENTS AND WASTE WATERS

( Fourth Reprint OCTOBER 1995 )

UDC 628.312:615.9:639.2

© Copyright 1972

**BUREAU OF INDIAN STANDARDS**  
MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG  
NEW DELHI 110002

## Indian Standard

# BIO-ASSAY METHODS FOR EVALUATING ACUTE TOXICITY OF INDUSTRIAL EFFLUENTS AND WASTE WATERS

Water Sectional Committee, CDC 26

### Chairman

DR T. R. BHASKARAN

### Members

PROF S. J. ARCEIVALA

DR P. R. BAVDEKAR

DR R. P. DAROGA ( Alternate )

SHRI A. K. BHATTACHARYA

SHRI S. DAS GUPTA ( Alternate )

CHEMIST AND METALLURGIST,  
CENTRAL RAILWAYS, BOMBAY

ASSISTANT RESEARCH OFFI-  
CER (CM-II), RDSO, CHIT-  
TARANJAN ( Alternate )

SHRI K. L. GHOSH

DR T. GUHA

DR R. S. VENKATARAMAN ( Alternate )

DR M. I. GURBAXANI

SHRI A. S. HARIHARAN

PROF B. L. MAHABAL

PROF N. MAJUMDAR

SHRI V. K. MALIK

SHRI K. KAMESHWARA RAO ( Alternate )

SHRI S. K. MITRA

MUNICIPAL ANALYST

DR S. D. PAUL

DR S. C. PILLAI

REPRESENTATIVE

SHRI T. VIRARAGHAVAN ( Alternate )

SHRI B. A. SHOLAPURWALA

SHRI V. M. SHAH ( Alternate )

SHRI PARITOSH C. TYAGI

SHRI M. R. VERMA

SHRI V. M. BHUCHER ( Alternate )

SHRI D. DAS GUPTA,

Director ( Chem )

### Representing

Geo-Miller & Co Private Ltd, Calcutta

Central Public Health Engineering Research Institute  
( CSIR ), Nagpur

Tata Chemicals Ltd, Bombay

National Test House, Calcutta

Railway Board ( Ministry of Railways )

Regional Research Laboratory ( CSIR ), Bhubaneswar  
Alkali & Chemical Corporation of India Ltd, Calcutta

The Tata Iron & Steel Co Ltd, Jamshedpur  
King Institute, Madras

Victoria Jubilee Technical Institute, Bombay

Indian Council of Medical Research, New Delhi

All India Distillers' Association, New Delhi

Ministry of Defence ( DGI )

Municipal Corporation of Greater Bombay, Bombay

Bhabha Atomic Research Centre, Bombay

Indian Institute of Science, Bangalore

Ministry of Health & Family Planning

Italab Private Ltd, Bombay

Planning Commission, New Delhi

National Physical Laboratory ( CSIR ), New Delhi

Director General, ISI ( Ex-officio Member )

### Secretary

DR G. M. SAXENA

Deputy Director ( Chem ), ISI

( Continued on page 2 )

**BUREAU OF INDIAN STANDARDS**

MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG

NEW DELHI 110002

(Continued from page 1)

**River Water and Industrial Effluents Subcommittee, CDC 26 : 1**

*Convener*

**PROF S. J. ARCEIVALA**

*Representing*

Central Public Health Engineering Research Institute  
( CSIR ), Nagpur

*Members*

**DR G. J. MOHANRAO** ( Alternate to  
Prof S. J. Arceivala )

**DR T. BHADURY**

Swasthya Sewa Nideshak, Government of Uttar  
Pradesh, Lucknow

**SHRI N. D. CHHABRIA**

Dorr Oliver India Ltd, Bombay

**SHRI R. NATARAJAN** ( Alternate )

**DIRECTOR ( CSM )**

Central Water & Power Commission, New Delhi

**SHRI N. C. RAWAL** ( Alternate )

**SHRI A. S. HARIHARAN**

King Institute, Madras

**SHRI D. V. KALRA**

Indian Chemical Manufacturers' Association, Calcutta

**SHRI V. K. DIKSHIT** ( Alternate )

**SHRI K. R. KRISHNASWAMI**

Madras Fertilizers Ltd, Madras

**SHRI B. C. MAJUMDAR**

Indian Leather Technologists' Association, Calcutta

**SHRI M. N. MOHANTY**

Orissa River Board and Director of Industries, Govern-  
ment of Orissa, Cuttack

**SHRI M. MAJEED MOHIUDDIN**

Indian Paper Mills Association, Calcutta

**SHRI B. N. SRIMANI**

Jadavpur University, Calcutta

**Panel for River Water and Waste Waters, CDC 26 : 1 : 11**

*Convener*

**DR T. R. BHASKARAN**

Geo-Miller & Co Private Ltd, Calcutta

*Members*

**DR T. BHADURY**

Swasthya Sewa Nideshak, Government of Uttar  
Pradesh, Lucknow

**DIRECTOR ( CSM )**

Central Water & Power Commission, New Delhi

**SHRI N. C. RAWAL** ( Alternate )

**SHRI T. DURAIRAJ**

Ministry of Health & Family Planning

**SHRI T. K. VEDARAMAN** ( Alternate )

**SHRI A. S. HARIHARAN**

King Institute, Madras

**PROF N. MAJUMDAR**

Indian Council of Medical Research, New Delhi

**SHRI M. MAJEED MOHIUDDIN**

Orient Paper Mills Ltd, Brajrajnagar

**DR A. PANDA** ( Alternate )

**SHRI K. C. PILLAI**

Bhabha Atomic Research Centre, Bombay

**DR N. U. RAO**

Central Public Health Engineering Research Institute  
( CSIR ), Nagpur

**SHRI R. P. MISRA** ( Alternate )

**SHRI G. S. RATNA**

Delhi Cloth & General Mills Co Ltd, New Delhi

**SHRI A. G. NENE** ( Alternate )

**SHRI M. A. SAMPATHKUMARAN**

All India Institute of Hygiene and Public Health,  
Calcutta

## *Indian Standard*

# BIO-ASSAY METHODS FOR EVALUATING ACUTE TOXICITY OF INDUSTRIAL EFFLUENTS AND WASTE WATERS

### 0. FOREWORD

**0.1** This Indian Standard was adopted by the Indian Standards Institution on 24 December 1971, after the draft finalized by the Water Sectional Committee had been approved by the Chemical Division Council.

**0.2** The bio-assay procedures given in the standard are intended for use in laboratories for evaluating the toxicity of industrial effluents and waste waters to fish. The test can be used to determine whether or not an effluent or effluent component is acutely toxic and, if toxic, the degree of toxicity. It also serves as a basis for judging whether or not an effluent can be discharged at a given rate without causing direct injury to fish and other aquatic organisms in the receiving water due to acute toxicity.

**0.3** The bio-assay is very useful in connection with the control of effluent disposal. Chemical examination of complex industrial effluents alone does not usually yield sufficient information. Many of their various toxic components cannot be readily detected, separated and measured by chemical means. Moreover, the degree of toxicity of each of the numerous substances and mixtures of chemicals is not known. The toxicity of effluents can be greatly influenced by interaction between their individual components and the dissolved minerals present in widely varying amounts in the receiving waters. Therefore, the toxicity of an industrial effluent to local fish in their natural habitat has to be evaluated directly through biological test under appropriate experimental conditions.

**0.4** In the preparation of this standard assistance has been taken from the following publications:

ASTM Designation: D 1345-59. Method of test for evaluation of acute toxicity of industrial waste water to fresh-water fish. American Society for Testing and Materials.

Standard methods for the examination of water and waste water. 1965. Ed 12. American Public Health Association, American Water Works Association, and Water Pollution Control Federation, New York.

## 1. SCOPE

**1.1** This standard lays down the bio-assay methods for evaluating the acute toxicity of industrial effluents and waste waters to fish.

## 2. GENERAL PRINCIPLES AND DIRECTIONS

**2.1 Measure of Acute Toxicity**—The measure of acute toxicity is the concentration of the test sample in a suitable diluent (referred to as diluent water) at which just 50 percent of the test animals are able to survive for a specified period of exposure. This is known as the median tolerance limit ( $TL_m$ ). The concentration of the sample in the diluent water is expressed as percent by volume. For example, a 10 percent concentration or a  $TL_m$  of 10 percent equals 1 part of the test sample in 9 parts of the diluent water. To evaluate toxicity, different concentrations of the test sample are tested so that the concentration lethal to 50 percent of the test animals (within the prescribed period) may be found or estimated by interpolation.

**2.2 Provision of Dissolved Oxygen During Test**—Since death of the fish may be caused also by deficiency of dissolved oxygen (D.O.) in the waste water, it is necessary that adequate dissolved oxygen is maintained in the sample under test.

**2.3 Uniformity of Experimental Procedure**—Reasonable uniformity of experimental procedure and of the manner of presentation of results is essential; widespread adoption of uniform methods will promote the accumulation of comparable data. However, rigid standardization of experimental material and conditions is not desirable, for it would tend to defeat the purpose of practical test. Strict comparability of test results could be achieved only by sacrificing much of their relevance and applicability to practical problems in specific localities. Experimental water (water used as a diluent) and test animals best suited for the purpose of each bio-assay should be selected. With respect to other features, the test can be more or less uniform.

**2.4 Sampling**—Samples of effluents or waste waters which are not constant in their composition shall be collected at different times and shall not be unnecessarily combined to make composite samples, because knowledge of the maximal toxicity of a variable effluent is often required in connection with the control of effluent disposal in flowing water rather than knowledge of the average toxicity. Extensive damage to the aquatic life of a receiving stream can result from brief, intermittent discharges of a highly toxic waste, even though the toxicity of the effluent at other times and its average toxicity are negligible. Therefore, if the composition and toxicity of an effluent vary considerably, it is necessary to test a number of individual (grab) samples taken at times when the effluent is likely to be most toxic, in order to determine the maximal toxicity. A composite sample of an effluent consisting of portions collected at regular time intervals can be useful only when average toxicity is to be evaluated.



**2.5 Storage of Samples** — Samples of industrial effluents or waste waters should be stored in completely filled, stoppered bottles at a fairly uniform temperature not greatly exceeding the initial temperature. Samples containing photosensitive material should be kept in brown bottles or in bottles covered with dark paper. If the effluent contains organic matter subject to bacterial decomposition, the samples shall be refrigerated (without freezing) and held at a temperature between 0°C and 4°C. Duration of storage prior to testing should be kept to a minimum.

**2.6 Interference** — Extreme volatility, instability and rapid detoxification of effluent components and excessive oxygen demand can render the Routine Bio-assay Method inapplicable, or can cause serious bias in the results. Extreme volatility, instability, or detoxification of important constituents is indicated when the recorded average survival time of the test animals in a fresh medium is much less than the survival time in the corresponding older, used medium, if adequate D.O. is present throughout both tests.

### 3. SELECTION OF METHOD OF TEST

**3.1** The basic Routine Bio-assay Method given under 4 constitutes the simplest procedure. It is widely applicable, being suitable for the detection and evaluation of acute toxicity which is not associated with excessive oxygen demand and which is due to substances that are relatively stable and are not extremely volatile. The Routine Bio-assay Method is designed so that the surface absorption of oxygen from the atmosphere plus some oxygen from the diluent generally provides an adequate amount of D.O. for the fish during the test period. Many industrial effluents have a high chemical or biochemical oxygen demand, which may cause oxygen depletion in test solutions. Although it is usually necessary to use artificial supplies of air or oxygen, D.O. shall nevertheless be maintained at levels adequate for the test fish.

**3.2** Uncontrolled aeration with compressed air has been generally unsatisfactory in bio-assays of industrial effluents. Small additions of air will not ordinarily maintain the necessary quantities of D.O., whereas vigorous aeration may drive off volatile toxicants and may greatly speed up biological oxidation.

**3.3** The following three methods of maintaining adequate D.O. are given under Modifications of Routine Bio-assay Method in 5:

- a) Controlled artificial oxygenation of test solutions,
- b) Initial oxygenation of diluent water, and
- c) Renewal of test solutions.

**3.4 Bio-assay Method Using *Daphnia*** — The microcrustacean *Daphnia magna* as an experimental animal in bio-assay may be used as an alternative procedure and has many advantages. These organisms are small in size (reaching a maximum size of 5 mm) and so a great number of them can be

reared in a small space. They can be grown in small bottles or cultured in mass in large aquaria. They are very sensitive to even low concentration of toxic substances. They are easy to culture and require only water containing bacteria or other equivalent for food. The organisms mature early giving birth to young ones within the first week of their life. After the first brood they give rise to a new brood every 2 or 3 days. An average of 20 to 30 may be produced in each brood. Their life span is short (about 2 months).

**3.5** When the experimental data are reported, any deviation from the basic routine procedure should be described in order to avoid misinterpretation of test results and their significance.

## **4. ROUTINE BIO-ASSAY METHOD**

### **4.1 Apparatus**

**4.1.1 Test Containers**—Six to twelve containers are usually required. These shall be chemically clean wide-mouth glass bottles or cylindrical glass jars 25 to 30 cm in diameter, capable of holding 15 to 20 litres of water. These dimensions would be suitable for tests with fish which are 5 to 7 cm long. Smaller or larger glass jars or rectangular glass aquarium tanks may be more suitable for smaller or larger fish.

**4.1.2 Acclimatizing Tanks**—Large rectangular glass aquarium tanks of 60 to 200 litres capacity, provided with arrangement for aeration with compressed or pumped air released near the bottom of the tank, shall be used. These shall be placed well in advance in the testing laboratory so that the temperature of the water may approximate the room temperature.

### **4.2 Test Animals**

**4.2.1 Types of Fish**—Fresh-water fishes commonly inhabiting unpolluted waters in the locality are preferable. They shall be adaptable to laboratory conditions of temperature, feeding and handling. The following fresh-water species are recommended:

- a) *Cyprinodontidae* (Top minnows)
  - 1) *Aplocheilichthys panchax*
- b) *Cyprinidae* (True minnows)
  - 1) *Nuria darinca*
  - 2) *Rasbora daniconius*
  - 3) *Danio Malabaricus*
  - 4) *Barbus ticto*
  - 5) *Barbus dorsalis*
  - 6) *Amblypharyngodon melettinus*
  - 7) *Cyprinus carpio*

Fish belonging to the same species shall be used for evaluation of each individual toxicity. They shall be identified as to genus and species.

**4.2.2 Source** — The fishes may be obtained from any single common source (stream, lake, hatchery, etc) but preferably from the body of water receiving the pollutant tested. They should all be collected and brought to the laboratory for acclimatization at about the same time.

**4.2.3 Size** — Specimens of more or less uniform size, averaging 5 to 7.5 cm in length are preferable. The largest shall be not more than 15 times as long as the smallest.

**4.2.4 Stocks** — Stocks of fish for test shall be kept initially in any suitable enclosures or containers (small ponds, concrete or wooden tanks, glass aquaria, etc) and in any water of suitable quality and temperature in which they will remain in good condition until needed.

**4.2.5 Acclimatization and Feeding** — The test fish shall be acclimatized for at least 10 days to laboratory conditions similar to those under which the tests are to be performed, with regard especially to the temperature and chemical characteristics of the dilution water. They shall be fed fairly regularly during the acclimatization period. They shall not be fed for about 48 hours before they are used in a test. They shall also not be fed during the test.

**4.2.6 Fitness** — The incidence of specimens dying or becoming seriously diseased in the acclimatizing aquarium during a period of 4 days immediately preceding a test shall be less than 10 percent. Otherwise the lot of test fish shall be deemed unfit for use. The specimens used for test shall not show any symptoms of disease or abnormalities of appearance or behaviour at the time of their transfer to the sample under investigation.

### **4.3 Diluent Water**

**4.3.1** The water to be used as a diluent and acclimatizing medium shall be obtained from the same body of water receiving the sample under investigation. When the toxicity of an effluent alone is to be determined the diluent water shall be obtained at a point where there is no pollution from any source.

**4.3.2** When the water receiving an effluent to be tested is subject to previous contamination with other effluents, the toxicity of the effluent under test in conjunction with the previous contamination shall be considered in judging rates of discharge which will not be acutely lethal. For this purpose, the diluent water shall be taken immediately upstream of the point of discharge of the effluent to be tested but outside the zone of its influence. Such an evaluation is possible only when the test animals can live in the diluent water.

**4.3.3** If uncontaminated diluent water cannot be obtained from the body of water under consideration, water of similar quality, with respect to its dissolved mineral content, should be obtained from another source,

or else prepared artificially from distilled or demineralised water to which appropriate amounts of minerals are added to simulate natural water conditions. The pH, temperature, alkalinity and hardness of this artificially prepared water shall be adjusted to match as closely as practicable those of the natural water.

**4.3.4** The diluent water shall be settled or filtered if excessive amount of suspended matter is present. The dissolved oxygen content of the diluent water shall be not below 4 mg/l. The diluent water shall be well aerated with compressed air. It shall not contain residual chlorine.

**4.4 Number of Test Fish** — At least 10 fish shall be used for each experimental concentration of the sample under test. They may be placed in one container with the test dilution, or preferably divided equally among two or more jars containing solutions of the same strength. It is desirable to maintain a ratio of about 1 g of weight of fish per litre of liquid.

**4.5 Dissolved Oxygen Content of Test Samples** — The D.O. content of test solutions, or dilutions, shall not fall below 4 mg/l throughout the period of the experiment.

#### **4.6 Procedure**

**4.6.1 Preliminary Test** — When testing effluents of completely unknown toxicity, much time and effort can be saved by conducting preliminary bio-assay to determine the approximate range of concentration of the effluent which should be covered in the full scale test. Prepare wide ranges of concentrations of the sample under test, for example, 100, 10, 1, 0.1 percent with the diluent water. Place two or more test fishes in an appropriate volume of each concentration contained in wide-mouth glass bottles. Observe after 24 hours. Select for full scale test the dilution ranges between the lowest concentration at which all fish survive for 24 hours and the highest concentration at which all or most fish die in 24 hours. The preliminary assay will also show whether excessive oxygen depletion occurs during the test period.

**4.6.2 Full Scale Bio-assay** — Test at least 4 to 6 concentrations arranged at equal logarithmic intervals within (and including) the limits of dilution ascertained from the preliminary test. A concurrent control test shall also be run under exactly similar conditions using the diluent water alone. The control shall be deemed satisfactory if not more than 10 percent of the fish placed in it die within the test period.

**4.6.3 Preparation of Test Dilutions of the Sample** — Shake the sample to disperse suspended matter uniformly. Withdraw measured portions of the sample and add to measured quantity of the diluent water and mix by gentle stirring. Avoid violent agitation.

**4.6.4 Transfer of Test Fish** — Transfer from the acclimatizing aquarium not less than 10 fish into each dilution. This shall be done within 30 minutes after preparation of the dilution. In transferring test fish special care shall be taken to handle them without causing injury.

#### **4.6.5 Duration of Test and Observations**

**4.6.5.1** The test fish shall not be fed during the course of the test.

**4.6.5.2** Observe the behaviour of the fish under test at regular intervals for at least 48 hours. The test period shall be extended to 96 hours if more than half of the test fish survive for 48 hours at the lowest concentration tested. Note especially the symptoms exhibited by the test fish such as abnormal behaviour and loss of equilibrium. The test fish shall be deemed dead if they show no movement even when gently prodded with a glass rod. The dead fish shall be removed as soon as observed.

**4.6.5.3** At the end of the test period transfer the surviving fish from each jar with as little agitation as possible to separate jars containing adequate amounts of the diluent water and observe for a further period of 24 hours. The number of fish which die in excess of those dying in the control within this period shall also be reported separately for the particular dilution.

**4.6.5.4** When required, minimum quantities of the test liquids may be removed for such determinations as dissolved oxygen. The volume of liquid in the test container shall then be made up with previously prepared dilution of similar strength.

**4.7 Estimation of  $TL_m$**  — Plot the percentages of surviving fish in each concentration on semi-logarithmic co-ordinate paper with the sample concentrations laid off on the logarithmic scale and survival percentages on the arithmetic scale. Fit a straight line to the points, particularly those between 20 and 80 percent survival. A line drawn from the point where this line crosses the 50 percent survival line to meet the concentration ordinate at right angles gives the  $TL_m$  concentration for the particular period of exposure involved.

#### **4.8 Reporting of Results**

**4.8.1** If all the test fish survive in the test sample (without any dilution) even after 96 hours the sample shall be reported as free from acute lethal toxicity.

**4.8.2** While reporting the  $TL_m$  value, mention shall be made of the species of fish used and their origin, weight and length; the temperature at which the test was carried out; pH and D.O. content of the test solution; the duration of the test; and the source of the diluent water and its mineral composition.

## 5. MODIFICATIONS OF ROUTINE BIO-ASSAY METHOD

**5.1 Controlled Artificial Oxygenation of Test Solutions**—If the sample under test has an excessive biochemical or other oxygen demand, it would be necessary to artificially aerate the test solution to maintain the minimum dissolved oxygen concentration. This shall be done by bubbling (30 to 180 bubbles per minute) pure oxygen from an oxygen cylinder provided with pressure reduction valve and arrangements for proper disbursement of the gas without violent agitation of the test animals. Supersaturation with dissolved oxygen shall be avoided.

**5.2 Initial Oxygenation of Diluent Water**—Another modification is the addition of oxygen to the diluent water. An industrial effluent can be free of oxygen and it can also have a considerable immediate (chemical) oxygen demand which shall be satisfied. The D.O. content of each test dilution shall be adequate at the beginning of the test. In order to introduce the required amount of oxygen without recourse to artificial aeration or excessive dilution of the test effluent, additional oxygen may be dissolved in the diluent water before it is mixed with the effluent. The oxygen content or the immediate oxygen demand of the effluent sample shall be determined. The required oxygen content of the diluent to be added, or the minimum degree of dilution with well-oxygenated (supersaturated) water of known oxygen content that will ensure an adequate initial D.O. level, can then be estimated. Oxygenation of the diluent water shall be done by bubbling compressed oxygen gas through it in a tall container.

**5.3 Renewal of Test Solution**—A third modification involves renewal of the liquids tested for the purpose of maintaining more or less uniform concentrations of any volatile and unstable toxic components and adequate D.O. content. This modification is recommended especially when there are reasons for believing that the toxicity of the liquids declines rapidly during the course of a test. Periodic renewal of the liquid tested, at daily or other convenient intervals necessary to maintain the test conditions, presents fewer difficulties. This can be accomplished by transferring the test animals quickly, by means of a dip net, to a test container with fresh liquid. Renewal of the test medium at intervals of 24 hours is often both convenient and sufficient, but renewal at shorter intervals, such as 12 or 8 hours, is sometimes necessary or advisable. Constant-flow systems, in which the diluted wastes or test solutions are renewed continually (that is, replaced with fresh mixtures flowing constantly into the experimental containers), may also be used.

## 6. METHOD USING DAPHNIA

### 6.1 Apparatus

**6.1.1 Culture and Test Containers**—These shall be made of glass and be chemically clean. Wide-mouth bottles of capacity 150 ml and of medium height may be used.

## 6.2 Test Animals

**6.2.1 Type of *Daphnia*** — The *Daphnia* shall be of species adaptable to the laboratory conditions of test and temperature, easily handled, available at all times of the year, present in normal streams and of value to fish either directly or indirectly. *Daphnia magna* meets all these requirements and hence is usually adopted as the test animal. In order to obtain reliable results the use of organisms 8 hours old or less is recommended.

**6.2.2 Number of Test Animals** — Fifty organisms shall be used for testing each concentration of the material under investigation.

**6.2.3 Rearing of Experimental Animals** — The experimental animals shall be reared in the laboratory in the culture medium. The medium shall be prepared as follows:

Mix 5 g of air-dried horse manure with 20 g of dried sandy muck in 1 litre of water. Allow the infusion to stand for two days and strain through bolting silk cloth. Allow the filtrate to stand for 4 to 6 days and use the clear supernatant.

**6.2.4 Culture of Test Animals** — Use 150-ml wide-mouth bottles. Place one mature female in each of a series of bottles filled with the culture medium. After 4 or 5 days add to each bottle, every alternate day, 1 mg of a suspension prepared by mixing 1 mg of dried yeast with 1 ml of water. Transfer the young ones from the bottles every day to a stock tank to prevent depletion of food materials to the mother. Occasionally add some yeast to the medium in the stock tank. Use a translucent viewing screen for making observations. A dark field illumination helps. Rotating the bottles while observing will ensure visibility of all animals in it.

**6.2.5 Preparation of Animals for Test** — Remove all young *Daphnias* from the culture bottles with a pipette 8 hours prior to the beginning of the test. Wash the *Daphnia* three times in diluent water (see 6.3) by transferring to a bottle containing the diluent water and allowing them to remain in the water for 5 minutes.

**6.3 Diluent Water** — Unpolluted water which is low in dissolved salts and hardness, with pH between 6.5 to 8.5 shall be used.

NOTE — *Daphnia* are very sensitive to chlorine and copper. These shall be absent in the diluent water.

## 6.4 Procedure

**6.4.1 Preliminary Test** — Prepare dilutions of the sample under test covering a wide range of concentrations, for example, 100, 10, 1, 0.1 and 0.01 percent with diluent water. Place 10 or more of the *Daphnia* in each concentration and observe for a period of 24 hours. From this determine the approximate concentration range required as in 4.6.1 and carry out the full scale test.

**6.4.2 Full Scale Bio-assay** — Set up six 150-ml wide-mouth bottles. Prepare the various concentrations of the sample under investigation using the diluent water. Concentrations in a series of dilutions using the logarithmic values are preferable. Keep the sixth bottle with the diluent water only as control. Transfer 50 *Daphnia* prepared for test as in 6.2.5 into each bottle. Observe and record the number of *Daphnia* which are immobilised in each test container at the end of 12 hours and 24 hours after their introduction and also at the end of 48 hours if the tests are continued beyond 24 hours.

**6.5 Estimation of  $TL_m$**  — From the observations obtained, obtain the value of  $TL_m$  concentration as prescribed in 4.7.

## **6.6 Reporting of Results**

**6.6.1** If most of the test animals survive in the test sample (without any dilution) even after 48 hours, the sample shall be reported as free from acute toxicity.

**6.6.2** While reporting the  $TL_m$  value, mention shall be made of the test animals used, the temperature at which the test was carried out, the duration of the test, and the source of diluent water and its mineral composition.



# BUREAU OF INDIAN STANDARDS

## Headquarters:

Manak Bhavan, 9 Bahadur Shah Zafar Marg, NEW DELHI-110002

Telephones : 331 01 31  
331 13 75

Telegrams : Manaksanstha  
(Common to all Offices)

## Regional Offices:

|   | Telephone              |
|---|------------------------|
| Central : Manak Bhavan, 9 Bahadur Shah Zafar Marg, NEW DELHI 110002       | 331 01 31<br>331 13 75 |
| *Eastern : 1/14 CIT Scheme VII M, V.I.P. Road, Maniktola, CALCUTTA 700054 | 37 86 62               |
| Northern : SCO 335-336, Sector 34-A, CHANDIGARH 160 022                   | 60 38 43               |
| Southern : C.I.T. Campus, IV Cross Road, MADRAS 600113                    | 235 23 15              |
| †Western : Manakalaya, E9 MIDC, Marol, Andheri (East), BOMBAY 400093      | 832 92 95              |

## Branch Offices:

|  |            |
|--|------------|
| 'Pushpak', Nurmohamed Shaikh Marg, Khanpur, AHMADABAD 380001                       | 30 13 48   |
| ‡Peenya Industrial Area, 1st Stage, Bangalore-Tumkur Road,,<br>BANGALORE 560058    | 39 49 55   |
| Gangotri Complex, 5th Floor, Bhadbhada Road, T.T. Nagar, BHOPAL 462003             | 55 40 21   |
| Plot No. 21 Satyanagar, BHUBANESHWAR 751007  | 40 36 27   |
| Kalaikathir Building, 6/48 Avanashi Road, COIMBATORE 641037                        | 21 01 41   |
| Plot No. 43, Sector 16 A, Mathura Road, FARIDABAD 121001                           | 8-28 88 01 |
| Savitri Complex, 116 G.T. Road, GHAZIABAD 201001                                   | 8-71 19 96 |
| 53/5 Ward No. 29, R.G. Barua Road, 5th By-lane, GUWAHATI 781003                    | 54 11 37   |
| 5-8-56C L.N. Gupta Marg, Nampally Station Road, HYDERABAD 500001                   | 20 10 83   |
| R 14 Yudhister Marg, C Scheme, JAIPUR 302005                                       | 38 13 74   |
| 117/418 B Sarvodaya Nagar, KANPUR 208005   | 21 68 76   |
| Seth Bhawan, 2nd Floor, Behind Leela Cinema, Naval Kishore Road,<br>LUCKNOW 226001 | 23 89 23   |
| Patliputra Industrial Estate, PATNA 800013   | 26 23 05   |
| C/o Smt. Sunita Mirakhur,,<br>66 D/C Annexe, Gandhi Nagar, JAMMU TAWI 180004       | —          |
| T.C. No. 14/1421, University P.O., Palayam, THIRUVANANTHAPURAM 695034              | 6 21 17    |

## Inspection Offices (With Sale Point):

|   |          |
|---|----------|
| Pushpanjali, 1st floor, 205-A, West High Court Road, Shankar Nagar Square,<br>NAGPUR 440010 | 525171   |
| Institution of Engineers (India) Building 1332 Shivaji Nagar,<br>PUNE 411005                | 32 36 35 |

\* Sales Office is at 5 Chowringhee Approach, P.O. Princep Street,  
CALCUTTA 700072 27 99 65

† Sales Office is at Novelty Chambers, Grant Road, BOMBAY 400007 309 65 28

‡ Sales Office is at 'F' Block, Unity Building, Narasimharaja Square,  
BANGALORE 560002 22 39 71